


# ACCESSORIES

FOR OPERATING TABLES



# ACCESSORIES

---

	<b>Pages</b>
<b>A</b> = Attachment clamps	03
<b>B</b> = General surgery	03 – 05
<b>C</b> = Urology / Gynaecology	05 – 06
<b>D</b> = Vascular surgery	06
<b>E</b> = DORO® products	07 – 08
<b>F</b> = Orthopaedics / Traumatology	09 – 11
<b>G</b> = Products for chick® 7300 & chick® LP	12 – 13
<b>H</b> = Paediatrics	14
<b>I</b> = Head rest, leg plates	14
<b>K</b> = Spine surgery	15
<b>L</b> = Positioning aids	16 – 17
<b>M</b> = Accessories for Stille imagiQ2	18

**7300 + arcus** = For chick® 7300 and schaeerer® arcus models

**LP + 7300** = For chick® LP and chick® 7300 models

**7300** = For chick® 7300 models

**LP** = For chick® LP models

**imagiQ2** = For Stille imagiQ2 models

## Accessories Schaeerer Medical for operating tables The solution for all surgical procedures

Schaeerer Medical accessories stand out due to their multiple applications as well as the first class quality of material and manufacturing. All metallic elements are manufactured of high quality stainless steel.

In addition to conventional accessories for operating tables we also offer you a complete line of gel pads and special positioning systems, which help prevent patients from injury or pressure related issues.



**90506 Radial setting clamp with fixing screw, rotating**  
Accepts any mounting posts 0.6"-0.7". Accessories can be adjusted to a variety of angles with the rotating clamp.



**73040 Rail clamp**  
Securely locks all accessories with 1" x 1/4" mounting blades anywhere on the side rails. Accessories can be mounted vertically or horizontally.



**74006 Socket XL**  
Accepts any round mounting posts up to 5/8"D and rectangular mounting posts up to 3/8" x 1". Clamp can rotate and lock for a variety of angles.



**90521 Armboard with rapid-action fixing clamp, with ball and socket joint**  
Armboard adjustable with one-hand operation through central locking lever. With fixing strap with Velcro fastener.



**74027 Socket**  
Attaches to all OR tables with standard 3/8"W x 1 1/8"H side rails. Fits round mounting posts up to 5/8"D.




**90621 Lateral armboard**  
Easy to use multi-positioning armboard. With its three articulating joints the arm can be positioned anywhere above or below the table.



**74036 Universal rail clamp**  
Accepts both 5/8"D and 1" x 3/8" flat bar mounting posts. Spring-loaded and attaches to OR table side rail.



**74028 Standard armboard**  
Without pad (requires 74016 2" armboard pad, Deluxe or 74025 2" armboard pad, Softcare). Snaps on to the rail. Trigger release and lock mechanism with 180° full range of motion. 6"W x 26.5"L.



**74053 Camlock Clamp**  
Attaches to rail even over multiple drapes. Accepts accessories with round posts up to 1/2"D and flat posts up to 1" x 1/4". Fully autoclavable.



**74031 Extra wide armboard**  
Without pad (requires 74024 2" extra wide armboard pad, Deluxe or 74026 2" extra wide armboard pad, Softcare). Snaps on to the rail. Trigger release and lock mechanism with 180° full range of motion. 11"W x 26.5"L.

**74030 Adjustable height extra wide armboard**  
Without pad (requires 74024 2" extra wide armboard pad, Deluxe or 74026 2" extra wide armboard pad, Softcare). Vertically adjusts to the height of the table top. Extra wide. Snaps on to the rail. Trigger release and lock mechanism with 180° full range of motion. 11"W x 26.5"L.

# ACCESSORIES



**70522 Cross arm support**  
Allows for arm positioning above the chest on supine procedures or the arm to be maintained at proper height for lateral procedures. Without socket.

**74006 Socket XL**



**74019 Deluxe arthroscopic leg holder**  
Offers optimal patient positioning with dual independent braces that freely rotate 360° and spread from 3"L to 11"W to fit any size leg, at any angle.

**74005 Disposable arthroscopic leg holder pads**  
Disposable premium pads for arthroscopic leg holder or deluxe arthroscopic leg holder (12/case).

**74006 Socket XL**



**74009 Multit-axis arm positioner**  
Versatile and comfortable arm positioning for a variety of surgical procedures. Ball socket locking system for every patient size and position.

**74006 Socket XL**



**74018 Nissen straps (pair)**  
Wrapped with a 1" thick pad to secure the patient regardless of the angle of the patient's position. The bariatric nissen straps, 72" long, fit most bariatric patients and have a weight capacity of 600 lbs.

**73040 Rail Clamp (purchase 2)**



**90511 Flexible Anaesthesia screen,**  
Provides a stable yet flexible platform for optimal anaesthesia screen positioning.

**90506 Radial setting clamp with fixing screw, rotating (purchase 2)**



**74037 Shoulder supports (1 pair)**  
The shoulder supports clamp quickly to any operating table and with its ball and socket adjustment mechanism conforms to any size patient.

**74006 Socket XL (purchase 2)**



**74040 Patient restraint strap 96" (2 hooks)**

**74045 Patient restraint strap 96" (no hooks)**

Stabilizes patient positioning for a variety of procedures. The patient restraint straps clamp quickly to standard surgical tables. Available with or without hooks.



**74048 Lateral Brace**

Without clamp. Simple and effective lateral patient positioners clamp quickly to any standard operating room table while the three swivel positioners conform to any size patient.

**74006 Socket XL (purchase 3)**



**74004 Arthroscopic leg holder**  
Spring-loaded braces slide into place and snap-lock. Clamshell design with Velcro strap that wraps completely over leg-holder and provides stable positioning for the leg.

**74005 Disposable arthroscopic leg holder pads**  
Disposable premium pads for arthroscopic leg holder or deluxe arthroscopic leg holder (12/case).


**74006 Socket XL**



**74032 Foot extension**

With pad. Without clamps. Can be used horizontally as a table extension and vertically as a foot support for reverse Trendelenburg positions to optimize surgical patient positioning.

**73040 Rail Clamp (purchase 2)**



**74044 Minor procedure table**  
For arm and hand surgery. Clamps and locks to standard U.S. operating room tables in seconds. Extremely durable, light weight, easy-to-use and radiolucent. Adjusts up or down to be level with the OR table with no leg attachment. 15"W x 30"L.

**B 21**



**74013 Stirrups (1 pair)**  
Trigger-grip release and lock mechanism provides complete control of lithotomy and abduction/adduction. Range of motion is - 35° low lithotomy to 90° high lithotomy. Adduction -9° and abduction 25°. 600 lbs patient weight capacity.

**C 26**

**74047 Stirrup clamp**  
Required accessory. Purchase two for use with Platinum Stirrups



**74011 Knee stabilizer radiolucent**  
Foot support easily mounts to operating room side rails. Letting the surgeon replace the entire knee while the patient's foot is held firmly in place.

**B 22**

**74012 Disposable foot pads (12/case)**

**74006 Socket XL**



**74014 Stainless steel transport trolley (Stirrups)**  
Optimizes OR storage and access to the platinum stirrups. Holds all required accessories and clamps.

**C 27**



**92006 X-ray cassette frame**  
Telescopic handle. Cassette frame with open inside. Telescopic slide-out: 30" x 45".

**B 23**

arcus



**74021 Sterile view leg drapes**  
Allows complete visibility of patient's legs and positioning during intra-operative adjustments in the lithotomy position. Latex free (10/case).

**C 28**




**90955 Intermediate plate**  
With US-side rails. Without cut-out.

**90956 Intermediate plate**  
With US-side rails. With cut-out.

**C 24**

7300



**71106 Uro drainage system**  
Surgical fluid containment system flexes to accommodate surgeon's access to patient during gynaecologic and urologic procedures. Includes mounting system and 5 Uro bags.

**71107 Uro bag, sterile**  
Sterile plastic Uro bag with drain hose - fits over Uro drainage system (10/case).

**71107NS Uro bag, non-sterile**  
Non-sterile plastic Uro bag with drain hose - fits over Uro drainage system (10/case).

**C 29**



**529-0629-00 Side rail extension**  
Attaches to any side rail to provide 18" additional space for attaching accessories.

**71077 Side rail extension set (1 pair)**

**C 25**



**92004 Light leg plate**  
For temporary support of patient's legs during patient transfer to the OR table and initial positioning prior to the procedure. Length: 33".

**C 30**

arcus

# ACCESSORIES



**73002 Articulating leg board for chick® 7300**  
For temporary support of patient's legs during patient transfer to the OR table and initial positioning prior to the procedure. Articulates from 0 to - 55°.

C 31 7300



**schaerer® arcus 90597 Carbon-fibre board 55 inches**,  
To be mounted on seat plate interface (instead of seat plate). Unobstructed radiolucency, max. patient weight: 396 lbs. With pad.

90615 Side rail

D 36 arcus



**73001 Pelvic Base**  
With pad. Short 8" table extension with side rails and perineal cut-out for use with the chick® 7300 for gyn/urological procedures or use with the 36" Carbon Fiber imaging board.

C 32 7300 + arcus



**chick® 7300 71130 Carbon fiber IA board, 36"**  
To be mounted to Pelvic Base. Unobstructed radiolucency, max patient weight 400 lbs. With pad.

**73004 Carbon fiber IA board 55"**  
To be mounted on standard interface. Unobstructed translucency, max patient weight 400 lbs. With pad.

D 37 7300



**71110 Foot switch**  
With connecting cable. For 3 table movements: Up / down, Trendelenburg / Reverse-Trendelenburg, lateral tilt left / right.

C 33 7300



**502-0204-00 7300 Hand control extension cord**  
Extends the reach of the chick® 7300 hand control by 5 feet.

D 38 7300



**92033 Foot switch**  
With connecting cable. For 3 table movements: Up / down, Trendelenburg / Reverse-Trendelenburg, lateral tilt left / right.

C 34 arcus



**10821 Removable side rail**  
Attachable 12" side rail for use with the chick® LP, Chick IOT and 36" and 55" Carbon Fiber Boards (71130, 73004).

D 39 LP



**schaerer® arcus 90596 Carbon fiber board, 33 inches**,  
To be mounted on seat plate interface (instead of seat plate). Unobstructed radiolucency, max. patient weight: 793 lbs. With pad.

**90612 Carbon fiber board, 33 inches**  
To be mounted on standard interface (instead of head or leg plate). Max patient weight 330 lbs. With pad.

90615 Side rail

D 35 arcus



**71129 Carbon-fibre armboard**  
With pad 2". For insertion between pad and carbon fibre board. Enables unobstructed X-ray imaging with arm outstretched from the heart to the fingers. Length pad: 27".

D 40 7300 + arcus



**1001.100-US DORO® QR3 headrest system, aluminium,**  
Set includes:  
**1001.001 QR3 Skull clamp**  
**3002-00 Swivel adaptor**  
**3001-00-US base unit**  
Features 3-pin rigid cranial fixation for prone, supine, lateral and sitting positions. Skull clamp is compatible with DORO® radiolucent base units and adaptors.

E 41



**3030-00 DORO® Headrest System, radiolucent**  
Features 3-point cranial fixation for cranial and cervical spine procedures using CT, DSA, X-ray and Fluoroscopy. Skull clamp, swivel and transitional members are completely radiolucent with no metal parts to interfere with imaging.  
Set includes:  
**3031-00 Base unit radiolucent**  
**3033-00 Swivel adaptor radiolucent**  
**3034-00 Skull clamp radiolucent**  
**3032-10 DORO® Long transitional member radiolucent (6")**  
**3032-20 DORO® Short transitional member radiolucent (5")**  
**3030-001 Storage case headrest system radiolucent**  
**3006-20 Disposable skull pins Titanium, Adult**

E 46



**1204.001 DORO® easy connect navigation adaptor, Stryker**  
Attaches directly to the navigation ready QRS Skull Clamp 1001.001

E 42



**3030-300 DORO® Intra-Op CT headrest system, radiolucent**  
Features extra bracing and transitional members to allow increased reach for optimal ISO center positioning in the scanner. Skull clamp, swivel adaptors, transitional members and bracing are completely radiolucent. Navigation tracking devices can be mount on the Quick-Rail of the skull clamp or on its starbursts. Compatible with most OR/surgical imaging tables and imaging systems.  
Set includes:  
**3031-00 Base unit radiolucent**  
**3033-10 Swivel adaptor radiolucent, short**  
**3034-00 Skull clamp radiolucent**  
**3032-10 DORO® radiolucent long transitional member (6") (2 pieces)**  
**3032-20 DORO® radiolucent short transitional member (5")**  
**3032-21 DORO® radiolucent transitional member, bracing**  
**3030-001 Storage case**  
**3006-20 Disposable skull pins Titanium, Adult**

E 47



**1204.002 DORO® easy connect navigation adaptor, BrainLab**  
Attaches directly to the navigation ready QRS Skull Clamp 1001.001

E 43



**1204.003 DORO® swivel horseshoe headrest for adults with adjusting lever**  
Universal non-invasive headrest designed for adults for prone or supine positioning. Interchangeable gel pads.

E 44




**3036-300 Intra-Op CT horseshoe RS headrest, RDL**  
Provides radiolucent cranial support in the prone or supine position for adults and children for CT imaging and Fluoroscopy. Connects direct to transitional member for chin-to-chest positioning. Permits vertical and lateral adjustments. Includes extension bar for skeletal traction.  
Set includes:  
**3036-70 Horseshoe RS headrest, RDL, adult w/ extension bar**  
**3031-00 Base unit, radiolucent**  
**3032-10 DORO® radiolucent long transitional member (6")**  
**3032-20 DORO® radiolucent short transitional member (5")**

E 48



**3008-00 DORO® swivel horseshoe pediatric headrest with adjusting lever**  
Universal non-invasive headrest designed for children for prone or supine positioning. Interchangeable gel pads.

E 45



no picture available for the moment

E 49

no picture available  
for the moment

E 50

**3033-64 DORO® navigation adaptor, Stryker, radiolucent**



E 55

**3005-00 DORO® reusable skull pins for adults SS (3 pieces).**

**3004-00 DORO® reusable skull pins for children SS (3 pieces).**

**3005-50 DORO® reusable skull pins for adults, titanium - MRI-safe (3 pieces).**

**3004-50 DORO® reusable skull pins for children, titanium - MRI-safe (3 pieces).**



E 51

**3033-65 DORO® navigation adaptor, Medtronic, radiolucent**



E 52

**3033-66 DORO® navigation adaptor, BrainLab, radiolucent**



E 53

**3033-50 DORO® interface skull clamp, radiolucent (connects BrainLab vario guide to RDL SC)**



E 54

**3006-00 DORO® disposable skull pins for adults (36 pieces - 12 procedures).**

**3006-10 DORO® disposable skull pins for children (36 pieces - 12 procedures).**

**3006-20 DORO® disposable skull pins for adults, titanium - MRI-safe (36 pieces - 12 procedures).**

**3006-30 DORO® disposable skull pins for children, titanium - MRI-safe (36 pieces - 12 procedures).**



**Orthopaedic- and traumatological accessory for surgery of the extremities**

schaerer® MIS-Extension is a flexible, easily operable system for the orthopaedic and traumatological surgery of the extremities. schaeerer® MIS-Extension combines maximum control and precise movements of the patient's leg with easy handling and pronounced flexibility accompanied by maximum stability. It is easy to be attached to the operating tables chick® 7300 or schaeerer® arcus 501 with help of the trolley.

**Basic set consisting of:**

- 90775 Trolley for schaeerer® MIS-Extension
- 90750 Hip rest simple left, with pad
- 90751 Hip rest simple right, with pad
- 90768 Hip rest simple, roll with pad
- 90753 Leg rest left, with pad
- 90754 Leg rest right, with pad
- 90755 Triangular plate left, with pad
- 90756 Triangular plate right, with pad
- 90777 Perineal rod with pad, MIS DAA
- 90779 Perineal rod with pad, diameter 200 mm for hip arthroscopy
- 90761 Seat plate, with pad
- 90758 Bar
- 90769 Extension shoe, white
- 90778 Stability support (only for chick® 7300)



**90757 Seat plate**  
With pad for schaeerer® MIS-Extension. The especially shaped seat plate made of carbon fibre material enables an unobstructed imaging with the X-ray image intensifier.



**Hip surgery direct anterior approach and hip arthroscopy**

**Necessary accessories:**

Art. no. 90757	1 piece
Art. no. 90758	2 pieces
Art. no. 90769	2 pieces
Art. no. 90777	1 piece
Art. no. 90753 or 90754	2 pieces
Art. no. 90755 or 90756	2 pieces

**Recommended accessories:**

Art. no. 90775	1 piece
Art. no. 90774 or 90776	1 piece

(only for chick® 7300)




**90758 Bar**  
For schaeerer® MIS-Extension. Cantilevered support of the bars (no floor support). Easy one-hand operation.



**90769 Extension shoe**  
White for schaeerer® MIS-Extension.

# ACCESSORIES



## 73018 Orthopedic Kit

Orthopedic fracture attachment provides everything needed for supine and lateral lower extremity fracture repairs. Carbon fiber spars, patented scalloped seat section and integrated gross traction rail allow for optimal imaging, hip area fall through and easy traction manipulation.

- 10823 Leg transfer plate with cushion
- 10210 Radiolucent perineal post
- 10362 Perineal post cushion (conductive)
- 10836 Offset adapter with 2" pad
- 10908 Well leg support (Black)
- 10913 Traction boot assembly (2 pieces)
- 73011 Traction unit extension (2 pieces)
- 71070 Side rail socket (2 pieces)
- 73007 Ortho pelvic base
- 529-0629-00 Side rail extension (1 piece)
- 73008 Orthopedic abductor bar (2 pieces)
- 529-0871-00 Ortho abductor bar clamp, long (1 piece)
- 73015 Articulating bracket
- 73016 Surgical drapery rod (2 pieces)
- 73023 Angled plate without cushion
- 73033 Ortho equipment cart

F 61

7300 + arcus



## 74010 Shoulder assist limb positioner

Designed to eliminate the need for OR staff to hold the arm during surgery. Maintains a low profile while delivering exceptional control and patient positioning. Mounts onto the side of any standard operating room table and provides a non-sterile base for which sterile components are inserted. The upper extremity positioner and control tube are steam autoclavable.

F 62



## 74043 Stainless steel transport dolly (limb positioner)

Storage system to hold all the control tubes and clamps and to vertically to reduce the storage footprint of the extremity positioners.

F 65



## 74042 Lateral shoulder holder

Shoulder positioner adjusts from 0-90° and both internal and external rotation can be easily controlled and locked in place. Traction settings from 0-15 lbs with simple knob adjustment. All components except support frame are autoclavable.

74053 Camlock Clamp (required accessory, purchase 1)

F 63



## 74017 Total knee replacement system 24"

Ball joint allows the surgeon to position the operative leg in extension, flexion, tilt and rotation. Locks in any position with a simple turn of the wrist.

74023 Total knee replacement system 30"

## 74020 Sterile TKR boot liner

Sterile boot liners for the total knee replacement system (12/case).

74053 Camlock Clamp (required accessory, purchase 1)

F 66



## 74041 Disposable hand/forearm wrap

Sterile hand/forearm wrap for use with shoulder assist limb or lateral shoulder holder (12/case).

F 64



F 67

**74051 Extremity positioner**  
For positioning the lower extremity for foot or ankle procedures. Sealed seams for easy clean-up. 16"L x 12"W x 7"H



F 71

**74022 Disposable face mask**  
Securely holds head during the procedure without interfering or compromising the patient's airway. Easily attaches to head frame with soft Velcro straps (12/case).



F 68

7300 + arcus

**74000 Integrated beach chair attachment for schaefer® arcus 601**


Integrated beach chair designed specifically for the schaefer® arcus 601. Chair articulation is managed directly from the table hand control. Double-ball joint headrest allows patient's head to tilt 30° in every direction. Shoulder pull-away panels provide excellent surgical exposure. 500 lbs patient weight capacity.

**74001 Powered beach chair for schaefer® arcus 501 and chick® 7300 - rail mounted**

Fits all surgical tables. Rail mounted with external motor and hand control for shoulder chair positioning. Double-ball joint headrest allows patient's head to tilt 30° in every direction. Shoulder pull-away panels provide excellent surgical exposure. 500 lbs patient weight capacity.


**74002 Easy lift beach chair for schaefer® arcus 501 and chick® 7300 - rail mounted**

Fits all surgical tables. Rail mounted with gas spring assist positioning mechanism. Double-ball joint head rest allows patient's head to tilt 30° in every direction. Shoulder pull-away panels provide excellent surgical exposure. 500 lbs patient weight capacity.



F 72

**13200 Montreal lateral positioning device**  
Modular lateral positioning device with stackable radiolucent cylinders to accommodate any size patient.



F 73

**13111 Disposable pad 6' roll**  
Polyfoam padding to wrap around radiolucent montreal lateral positioning posts.

**13115 Disposable pad 6' roll (3 rolls)**



F 69

**74007 Stainless steel transport trolley (beach chair)**  
The durable, easy-to-use, transport dolly is designed to store any shoulder chair from Schaefer Medical. Stainless steel handle and frame as well as four caster wheels permit efficient movement and storage. Wide base resists tipping.



F 74

**10522 Disposable pad set with 24 x 9" perineal post pads**  
Complete disposable foam pad set for lower extremity fracture repair. Includes planter plate pad (24 pieces), heel rest pad (24 pieces), 9" perineal post pad (24 pieces).



F 70

**74003 Disposable head and chin straps**  
Holds head securely during the procedure. Velcro straps are soft and pliable and attach to head frame (12/case).



F 75

**10347 Disposable pad set with 12 saddle pads and 12 x 7" perineal post pads**  
Deluxe complete foam pad set. Includes saddle pad (12 pieces), 7" perineal pad (12 pieces), planter plate pad (24 pieces), heel rest pad (24 pieces) – All are foam (12/case).

# PRODUCTS FOR CHICK® 7300 + CHICK® LP



**10362 Perineal post pad - rubber**  
Replacement post pad – firm rubber.

G 76 LP + 7300



**10940 3/4" Disposable traction boot insert**  
¾ foam traction boot inserts (5/case).

G 81 LP + 7300



**10373 7" Disposable perineal post pad**  
Disposable foam 7" perineal post pads (12/case).

G 77 LP + 7300

**10374 9" Disposable perineal post pad**  
Disposable foam 9" perineal post pads (12/case).



**10210 Vertical perineal post without pad**  
Radiolucent perineal post. Can be used with all perineal post pads.

G 82 LP + 7300

**10828 Pediatric perineal post**



**10378 Disposable planter plate pad**  
(12 pairs/case)

G 78 LP + 7300



**10836 Offset adaptor**  
For use with chick® LP or 7300 Ortho kit to provide greater hip exposure for small patients and reduced direct perineal pressure.

G 83 LP + 7300



**10396 Lottes ankle pad insert**  
Foam ankle pad inserts (12 pairs/case).

G 79 LP + 7300



**10822 GM support board**  
For lateral procedures on the chick® LP. With pad.

G 84 LP



**10607 chick® LP isolation drapes with 3 pouches**  
Sterile isolation drapes with 3 pouches specifically designed for use with the chick® LP fracture table (12/case).

G 80 LP



**10824 Counter traction horizontal set**  
For lateral procedures on the chick® LP.

G 85 LP



**10531 Bell foot piece with traction**

**G 86** LP + 7300



**22111 90 Degree pin and wire set with ferrules**  
For use with chick® LP of chick® 7300 ortho kit for skeletal traction. Includes ferrule set. Attaches to fine traction barrel.

**G 91** LP + 7300



**10825 Tibia extender**  
Provides knee flexion for tibia repair.

**G 87** LP



**Chick roller**  
Ergonomic method for staff to transfer patients to and from stretchers, surgical tables and hospital beds with minimal effort and maximum safety. Heavy duty 5 cylinder construction with static prevention cover. Available in two sizes.

**13001 Standard roller with conductive cover**  
1 ¼" x 14 ¼" x 67"

**13030 Short roller with conductive cover**  
1 ¼" x 14 ¼" x 30"

**13033 Wall mounting kit for Chick Roller**

**G 92** LP + 7300



**10908 Articulating well leg support**  
Well leg support system that adjusts in all directions to fit all patient sizes. Ball and joint adjustment of the well leg support pad for exact positioning.

**G 88** LP + 7300



**10005 Infant cast table**  
Designed specifically for use with children. It can be used for full body casting, hip or shoulders SPICAS and can be placed on a surgical table, stretcher, exam table or crib.

**G 93** LP + 7300



**10913 Traction boot assembly**  
For use with the chick® LP or chick® 7300 ortho kit. Includes fine traction control barrel and boot assembly with ratchet tightening mechanism.

**G 89** LP + 7300



**10006 Infant cast table with traction**  
Includes the infant cast table and removable traction system.

**G 94** LP + 7300



**22313 Traction cuff assembly**  
For use with the chick® LP or chick® 7300 ortho kit. Includes fine traction control barrel and boot assembly with air bladder for foot capture.

**G 90** LP + 7300

# ACCESSORIES



**90617 Child surgery plate**  
With US-side rails. Child surgery plate for using on models chick®7300 and schaefer® arcus 501/701. With pad. Max. patient weight 77 lbs.

H 95 7300 + arcus



**92065 Leg plate, one piece, tilting up and down**  
With pad. Easily articulates 15° to -90° and is infinitely positionable with gas spring lift assist.

I 100 arcus



**92410 Head rest, standard, tilting up and down**  
Unlimited easy positioning from 40° to -40° with gas spring assist.

I 96 arcus




**74033 Table width extension 20"L x 4"W**  
Bariatric table width extension mounts on the OR table side rail. Unique design allows any table accessory to be mounted to this extension.

I 101



**92411 Head rest, double joint, tilting up and down**  
The headrest can be manipulated in two directions. Separate gas spring lift assist for the whole assembly and the adjustable top.

I 97 arcus



**74034 Table width extension 15"L x 8"W**  
Bariatric table width extension mounts on the OR table side rail. Unique design allows any table accessory to be mounted to this extension.

I 102



**92055 Leg plates, 2-part, tilting up and down**  
Split legs consist of two small and two large leg plates which can be moved in multiple directions. With gas spring lift assist.

I 98 arcus



**74035 Table width extension 11"L x 8"W**  
Bariatric table width extension mounts on the OR table side rail. Unique design allows any table accessory to be mounted to this extension.

I 103



**92063 Leg plate, one piece**  
Without articulation, with pad, with US-side rails 28" x 34" x 4".

I 99 arcus



**92431 Intermediate plate, long**  
With US-side rails.

**92432 Intermediate plate, short**  
With US-side rails.

I 104 arcus



**74055 SchureSpine Frame**  
 The radiolucent Spine Frame provides a stable and convenient platform for spine surgeries. It can be used with any O.R. table. Adjustable for optimum exposure. Patient weight capacity of 500 lbs (227 kgs).

- Includes:
- Tapering adjustable softcare pads.
  - Arm Cradle, disposable (pair)
  - Slotted Head Rests, disposable (each)
  - Crank Handle

**K 105**

## OPTIONAL ACCESSORIES



**74049 Disposable Slotted Arm Cradles**  
 (12/case), 5"Wx24"Lx3"H

**K 106**



**74059 Spine Frame Spine Dolly**

**K 108**



**74050 Disposable Slotted Head Rest**  
 (6/case), 8 1/2"L x 8"W x 3 3/4"H  
 (22cm x 20cm x 10cm)

**K 107**

**76003 Spine Frame Disposable Kit**  
 Includes 2 paper pad covers 2 foam bar covers  
 (20 Set/case).

# POSITIONING AIDS



## Patient positioning with vascocare Evolution®

The unique patient positioning system:

- Relieves pressure and stabilizes through anatomically formed viscoelastic elements - also in Trendelenburg and tilt positions.
- Minimum influence on blood circulation through optimum distribution of pressure as evidenced with "FLP" blood pressure measurements.
- Minimum influence on skin temperature and moisture through optimized materials.
- Low skin friction / shear forces through sliding, skin friendly surfaces.
- Robust and simple handling with wear resistant, wipe-disinfectable surfaces.



**75002 Head support**  
Supine position, version 1 - small (pediatrics).

**75003 Head support**  
Supine position, version 1 - large.

L 109



**75004 Head support**  
Supine position, version 2 - medium.

**75005 Head support**  
Supine position, version 2 - large.

L 110

### Disposable covers for supine head positioners:

**75006 Disposable covers for 75002,**  
25 pieces.

**75007 Disposable covers for 75003,**  
25 pieces

**75008 Disposable covers for 75004,**  
25 pieces

**75009 Disposable covers for 75005,**  
25 pieces

**75010 Disposable Measuring Tape,**  
50 pieces

**XXXXX Strap and Retainer Kit**



L 111

### 521300 Torso positioners

Supine position, large (1 pair)

**XXXXXX Torso positioner**

Supine position, large (left)

**XXXXXX Torso positioner**

Supine position, large (right)



L 112

### 521400 Knee positioners

Supine position, large (1 pair)



L 113

### 521500 Foot positioners

Supine position, large (1 pair)

**XXXXXX Foot positioner**

Supine position, large (left)

**XXXXXX Foot positioner**

Supine position, large (right)



L 114

### 75000 vascocare Evolution® SET

Head support, version 2 – large

Torso positioners, supine position (1 pair)

Knee positioners, supine position (1 pair)

Foot positioners, supine position (1 pair)

### 75001 vascocare Evolution® deluxe SET

Head support, version 2 – large

Torso positioners, supine position (1 pair)

Knee positioners, supine position (1 pair)

Foot positioners, supine position (1 pair)

Arm positioners, supine position (1 pair)

**GEL-Pads**

All Schaerer® GEL pads are filled with a special POWERGEL, offering optimum pressure distribution. The special surface is wear resistant, gentle to the skin and has sliding properties. The Consistens-Gel is not flowing. It only adapts its form.



**L 115**

**Head support, donut head pads**

**510001** Dimensions: 4" x 2" x 1"


**510003** Dimensions: 8" x 3" x 2"



**L 116**

**Head support, horseshoe head pads**

**510014** Dimensions: 8" x 3" x 3"



**L 117**

**Table pads, in one piece**

**510302** Dimensions: 39" x 20" x 0.4"

**510403** Dimensions: 46" x 20" x 0.4"



**L 118**

**74039 Double angle wedge 24"x9-3/4"x 6-3/4"**

Positioning pads are reusable and ideal for surgical extremity positioning. Product is radiolucent and is coated with argon antimicrobial ionic silver.



**L 119**

**74049 Disposable slotted arm cradles**

Disposable high-density polyfoam arm cradles are supportive, comfortable and economical (12/case).



**L 120**

**74050 Disposable slotted head rest**

Disposable, high-density polyfoam head rests are supportive, comfortable and economical (6/case).

# ACCESSORIES FOR STILLE IMAGIQ2 PRODUCTS



**535-1720 Side rails US**  
Attaches directly to carbon fibre table top to allow for additional rail mounted accessories.



**535-1760 Pan handle**  
Locks, unlocks and controls the XY movement of the carbon fiber table top.



**535-1730 Fistula armboard**  
Detachable radiolucent armboard which allows for unobstructed imaging from the chest to the arm. No metal in the imaging area. Attaches directly to the carbon fibre table top.



**535-1761 Foot control**  
Allows for foot control adjustment of table height, iso-centric roll, and Trendelenburg / Reverse Trendelenburg.



**535-1740 Head rest adjustable**  
Anesthesia friendly head support is detachable and adjusts for convenient intubations and patient access.



**535-1741 Head rest fixed**  
Rigid head rest attachment for better access to the patient's head and elongation of the body imaging area.



**535-1751 Catheter tray 59"**  
**535-1750 Catheter tray 39"**  
**535-1752 Catheter tray 20"**  
Stainless steel catheter tray extension adds a functional work surface for guide wires, instruments and other tools.



Available in California, Nevada and Hawaii through Covenant  
Call: 1(800) 975-0563  
Email: [Info@covenant-hcec.com](mailto:Info@covenant-hcec.com)  
Visit: [www.covenant-hcec.com.com](http://www.covenant-hcec.com.com)